

SPECIFICATION

Perillaldehyde

for the Japanese Pharmacopoeia Crude Drugs Test (for Assay and TLC)

REQUIREMENT	SPECIFICATION
Description (JP method)	Colorless to light brown transparent liquid
Absorbance E1%1cm (230nm) (JP method)	850~950
Identification (JP method)	to pass test
Purity (JP method)	-
Related substances (1) (JP method)	to pass test
Related substances (2) (JP method)	to pass test
* Solubility in methanol	to pass test
* Water	max.1.0%

This product is suited for Perillaldehyde for assay and thin-layer chromatography of general tests on The Japanese Pharmacopoeia. (*:Additional test performed by Wako)