

## SPECIFICATION

## Potassium Dichromate

TraceSure(R)

REQUIREMENT	SPECIFICATION
* Appearance Assay (K <sub>2</sub> Cr <sub>2</sub> O <sub>7</sub> ) Water-insoluble matter Chloride (Cl) Sulfate (SO <sub>4</sub> ) Sodium (Na) Copper (Cu) Calcium (Ca) Lead (Pb) Iron (Fe)	Reddish yellow ~ yellowish red, crystals ~ crystalline powder min.99.98%(mass/mass)  max.0.005%(mass/mass)  max.5mg/kg  max.0.03g/kg  max.0.1g/kg  max.1mg/kg  max.0.02g/kg  max.1mg/kg  max.2mg/kg

JIS K 8005:2014\*:Additional test performed by Wako but not required by JIS