

## SPECIFICATION

3 PFCs Mixture Standard Solution (PFHxS,PFOS,PFOA each 2  $\mu\text{g}/\text{mL}$  Methanol Solution)  
for Water Analysis

REQUIREMENT	SPECIFICATION
Appearance	Colorless clear liquid
Practical test	to pass test
Perfluorooctanoic acid (PFOA) (linear chain)	1.97~2.03 $\mu\text{g}/\text{mL}$
Perfluorohexanesulfonic acid (PFHxS) (linear chain)	1.97~2.03 $\mu\text{g}/\text{mL}$
Perfluorooctanesulfonic acid (PFOS) (linear chain)	1.97~2.03 $\mu\text{g}/\text{mL}$