

Code No.19-0010

SPECIFICATION

Salicylaldehyde

Guaranteed Reagent

REQUIREMENT	SPECIFICATION
* Appearance	Colorless ~ pale yellow, clear liquid
Assay (C ₇ H ₆ O ₂) (GC)	min.97.0%(mass/mass)
Solubility in ethanol	to pass test
Density (20°C)	1.163~1.168g/ml
Refractive index 20°C	1.571~1.574
Water	max.0.2%(mass/mass)
Suitability for acetone determination	to pass test

JIS K 8390:2011

(*:Additional test performed by Wako but not required by JIS)