

Code No.19-0106

SPECIFICATION

Sodium Acetate Trihydrate

Wako 1st Grade

REQUIREMENT	SPECIFICATION
Appearance	White, crystals ~ crystalline powder
Solubility in water	to pass test
pH of a 50g/L solution at 25°C	7.5~9.0
Chloride (Cl)	max.0.02%
Sulfate (SO ₄)	max.0.01%
Copper (Cu)	max.0.001%
Lead (Pb)	max.0.001%
Iron (Fe)	max.0.002%
Substances reducing permanganate (as O)	max.0.004%
Assay	min.98.5%