

SPECIFICATION

Sodium Acetate

Guaranteed Reagent

REQUIREMENT	SPECIFICATION
* Appearance	White, crystals ~ crystalline powder
Assay (CH ₃ COONa)	min.98.5%(mass/mass)
Solubility in water	to pass test
pH of a 50g/l solution at 25°C	7.5~9.0
Chloride (Cl)	max.0.005%(mass/mass)
Phosphate (PO ₄)	max.0.001%(mass/mass)
Sulfate (SO ₄)	max.0.003%(mass/mass)
Potassium (K)	max.0.02%(mass/mass)
Copper (Cu)	max.5ppm(mass/mass)
Calcium (Ca)	max.0.005%(mass/mass)
Lead (Pb)	max.0.001%(mass/mass)
Iron (Fe)	max.0.001%(mass/mass)
Substances reducing permanganate	max.0.002%(mass/mass)

JIS K 8372:2013

(*:Additional test performed by Wako but not required by JIS)