

Code No.19-0174

# SPECIFICATION

## Sodium Chromate Tetrahydrate

Wako Special Grade

REQUIREMENT	SPECIFICATION
Appearance	Yellow, crystals ~ crystalline powder
Solubility in water	to pass test
pH of a 10g/L solution at 25°C	8.0~10.0
Chloride (Cl)	max.0.005%
Sulfate (SO <sub>4</sub> )	max.0.01%
Potassium (K)	max.0.05%
Copper (Cu)	max.0.001%
Calcium (Ca)	max.0.02%
Lead (Pb)	max.0.001%
Iron (Fe)	max.0.002%
Assay	min.99.0%