

Code No.19-0179

SPECIFICATION

Trisodium Citrate Dihydrate

Wako 1st Grade

REQUIREMENT	SPECIFICATION
Appearance	White, crystals ~ crystalline powder
Solubility in water	to pass test
pH of a 50g/L solution at 25°C	7.5~9.0
Chloride (Cl)	max.0.01%
Sulfate (SO ₄)	max.0.02%
Nitrogen compounds (as N)	max.0.005%
Copper (Cu)	max.5 ppm
Calcium (Ca)	max.0.02%
Lead (Pb)	max.0.001%
Iron (Fe)	max.0.001%
Assay	min.99.0%