

## SPECIFICATION

## Sodium Formate

## Guaranteed Reagent

REQUIREMENT	SPECIFICATION
* Appearance Assay (HCOONa) Solubility in water pH of a 50g/l solution at 25°C Chloride (Cl) Phosphate (PO <sub>4</sub> ) Sulfate (SO <sub>4</sub> ) Magnesium (Mg) Calcium (Ca) Lead (Pb) Iron (Fe)	White, crystals ~ crystalline powder min.98.0%(mass/mass) to pass test 7.0~8.5 max.0.001%(mass/mass) max.0.001%(mass/mass) max.0.001%(mass/mass) max.0.001%(mass/mass) max.0.001%(mass/mass) max.5ppm(mass/mass) max.5ppm(mass/mass)

JIS K 8267-2014

(\*:Additional test performed by Wako but not required by JIS)