

SPECIFICATION

Sodium Dihydrogenphosphate Dihydrate

Guaranteed Reagent

REQUIREMENT	SPECIFICATION
* Appearance	White, crystals ~ crystalline powder
Assay (NaH ₂ PO ₄ ·2H ₂ O)	99.0~102.0%(mass/mass)
Solubility in water	to pass test
pH of a 50g/l solution at 25°C	4.1~4.5
Chloride (Cl)	max.0.001%(mass/mass)
Nitrate	to pass test
Sulfate (SO ₄)	max.0.005%(mass/mass)
Heavy metals (as Pb)	max.0.001%(mass/mass)
Calcium (Ca)	max.0.02%(mass/mass)
Arsenic (As)	max.1ppm(mass/mass)
Iron (Fe)	max.0.001%(mass/mass)
Ammonium (NH ₄)	max.0.002%(mass/mass)

JIS K 9009:2012

(*:Additional test performed by Wako but not required by JIS)