

## SPECIFICATION

## Sodium Diphosphate Decahydrate

## Guaranteed Reagent

REQUIREMENT	SPECIFICATION
* Appearance	White, crystals ~ crystalline powder
Assay (Na <sub>4</sub> P <sub>2</sub> O <sub>7</sub> ·10H <sub>2</sub> O)	min.99.0%(mass/mass)
Solubility in water	to pass test
pH of a 50g/l solution at 25°C	10.0~10.7
Chloride (Cl)	max.0.001%(mass/mass)
Nitrate	to pass test
Orthophosphate	to pass test
Sulfate (SO <sub>4</sub> )	max.0.005%(mass/mass)
Heavy metals (as Pb)	max.5ppm(mass/mass)
Arsenic (As)	max.1ppm(mass/mass)
Iron (Fe)	max.0.001%(mass/mass)
Ammonium (NH <sub>4</sub> )	max.0.001%(mass/mass)

JIS K 8785:2012(\*:Additional test performed by Wako but not required by JIS)