

## SPECIFICATION

## Sodium (+)-Tartrate Dihydrate

## Guaranteed Reagent

REQUIREMENT	SPECIFICATION
* Appearance Assay(C <sub>4</sub> H <sub>4</sub> Na <sub>2</sub> O <sub>6</sub> ·2H <sub>2</sub> O) Solubility in water pH of a 50g/L solution at 25°C Chloride (Cl) Phosphate (PO <sub>4</sub> ) Sulfur compounds (as SO <sub>4</sub> ) Copper (Cu) Calcium (Ca) Zinc (Zn) Lead (Pb) Iron (Fe) Ammonium (NH <sub>4</sub> )	White, crystals ~ crystalline powder min.99.0%(mass/mass) to pass test 7.0~9.0 max.0.002%(mass/mass) max.0.002%(mass/mass) max.0.005%(mass/mass) max.1ppm(mass/mass) max.0.005%(mass/mass) max.2ppm(mass/mass) max.1ppm(mass/mass) max.0.001%(mass/mass) max.0.002%(mass/mass)

JIS K 8540:2016

(\*:Additional test performed by Wako but not required by JIS)