

Code No.19-0413

## SPECIFICATION

## Strontium Acetate Hemihydrate

Wako Special Grade

REQUIREMENT	SPECIFICATION
Appearance	White crystalline powder
Solubility in water	to pass test
pH of a 50g/l solution at 25°C	6.0~9.0
Chloride (Cl)	max.0.001%
Nitrate	to pass test
Sulfate (SO <sub>4</sub> )	(NO <sub>3</sub> :max.about 0.005%) max.0.001%
Heavy metals (as Pb)	max.5ppm
Sodium (Na)	max.0.01%
Potassium (K)	max.0.01%
Calcium (Ca)	max.0.02%
Barium (Ba)	max.0.02%
Iron (Fe)	max.5ppm
Assay	min.99.0%