

SPECIFICATION

Amidosulfuric Acid

Guaranteed Reagent

| REQUIREMENT | SPECIFICATION |
|--|--------------------------------------|
| * Appearance | White, crystals ~ crystalline powder |
| Assay (HOSO ₂ NH ₂) | min.99.5%(mass/mass) |
| Solubility in water | to pass test |
| Residue after ignition (as sulfate) | max.0.005%(mass/mass) |
| Chloride (Cl) | max.5ppm(mass/mass) |
| Copper (Cu) | max.5ppm(mass/mass) |
| Lead (Pb) | max.5ppm(mass/mass) |
| Iron (Fe) | max.5ppm(mass/mass) |

JIS K 8587 : 2018

*:Additional test performed by Wako but not required by JIS