

## SPECIFICATION

Sulfuric Acid  
Guaranteed Reagent

REQUIREMENT	SPECIFICATION
* Appearance Assay (H <sub>2</sub> SO <sub>4</sub> ) Residue after ignition (as sulfate) Chloride (Cl) Nitrate Copper (Cu) Lead (Pb) Arsenic (As) Selenium (Se) Iron (Fe) Ammonium (NH <sub>4</sub> ) Substances reducing permanganate	Colorless clear liquid min.95.0%(mass/mass) max.5ppm(mass/mass) max.0.2ppm(mass/mass) to pass test max.0.1 ppm(mass/mass) max.0.1 ppm(mass/mass) max.0.02ppm(mass/mass) max.0.5ppm(mass/mass) max.0.5ppm(mass/mass) max.1 ppm(mass/mass) to pass test

JIS K 8951:2006

(\*:Additional test performed by Wako but not required by JIS)