

Code No.19-0553

## SPECIFICATION

### Trisodium Citrate Dihydrate

for Amino Acid Automated Analysis

REQUIREMENT	SPECIFICATION
Appearance	White, crystals ~ crystalline powder
Solubility in water	to pass test
Water-insoluble matter	max.0.001%
pH of a 40g/l solution at 25°C	7.8~8.6
Nessler reactive substances (as NH <sub>4</sub> )	max.0.3ppm
Ninhydrin reactive substances	to pass test
Assay	min.99.0%