

Code No.19-0576

## SPECIFICATION

### Sodium Isoascorbate Monohydrate

Wako Special Grade

REQUIREMENT	SPECIFICATION
Appearance	White ~ slightly yellow, crystalline powder
Solubility in water	to pass test
Specific rotation [ $\alpha$ ]D20(c=2,H2O)	+95~+98°
pH of a 20g/l solution at 25 °C	6.5~8.0
Heavy metals (as Pb)	max.0.002%
Arsenic (As)	max.2ppm
Assay	min.98.0%