

SPECIFICATION

5-Sulfosalicylic Acid Dihydrate
for Biochemistry

REQUIREMENT	SPECIFICATION
Appearance	White, crystals ~ crystalline powder
Solubility in water	to pass test
Solubility in acetic acid	to pass test
Heavy metals (as Pb)	max.5ppm
Iron (Fe)	max.5ppm
Salicylic acid (HOC ₆ H ₄ COOH)	max.0.01%
Suitability for cholesterol determination	to pass test
Assay	min.99.0%