

Code No.19-0749

# SPECIFICATION

## Disodium Hydrogenphosphate 12-Water for Oxidant Analysis for Oxidant Analysis

| REQUIREMENT                           | SPECIFICATION  |
|---------------------------------------|----------------|
| * Appearance                          | White crystals |
| Assay                                 | 99.0~100.5%    |
| Solubility in water                   | to pass test   |
| pH of a 50g/L solution at 25°C        | 9.0~9.4        |
| Nitrogen compounds (as N)             | max.0.001%     |
| Heavy metals (as Pb)                  | max.5ppm       |
| Suitability for oxidant determination | to pass test   |

JIS K 8828:2004

(\*:Additional test performed by Wako but not required by JIS)