

Code No.19-0866

SPECIFICATION

Saccharin

Wako Special Grade

REQUIREMENT	SPECIFICATION
Appearance	White ~ nearly white, powder
Solubility in water	to pass test
Solubility in ethanol	to pass test
Loss on drying at 110°C	max.1.0%
Chloride (Cl)	max.0.01%
Sulfate (SO ₄)	max.0.01%
Heavy metals (as Pb)	max.0.002%
Benzoic acid and salicylic acid	to pass test
Substances darkened by sulfuric acid	to pass test
Assay	min.98.0%