

## SPECIFICATION

## Stearic Acid

## Guaranteed Reagent

REQUIREMENT	SPECIFICATION
* Appearance Assay (C <sub>18</sub> H <sub>36</sub> O <sub>2</sub> ) (GC) Melting point Residue after ignition (as sulfate) Acid value Iodine value	White shot min.95.0%(mass/mass) 68~71°C max.0.01%(mass/mass) 195~200 max.1.5

JIS K 8585:2015(\*:Additional test performed by Wako but not required by JIS)