

Code No.19-1109

SPECIFICATION

Sodium Azide
Guaranteed Reagent

REQUIREMENT	SPECIFICATION
* Appearance	White ~ nearly white, crystalline powder or mass min.98.0%(mass/mass) to pass test max.0.2%(mass/mass) max.0.1%(mass/mass) max.2ppm(mass/mass)
Assay (NaN ₃)	
Solubility in water	
Alkalinity (as NaOH)	
Loss on drying at 105°C	
Chromium (Cr)	

JIS K 9501 : 2019(*:Additional test performed by Wako but not required by JIS)