

SPECIFICATION

Schizandrin

for the Japanese Pharmacopoeia Crude Drugs Test (for Thin-Layer Chromatography)

REQUIREMENT	SPECIFICATION
Description (JP method)	White crystals or crystalline powder
Melting point (JP method)	130~135°C
Purity (JP method)	-
Related substances (JP method)	to pass test
* Solubility in methanol	to pass test
* Loss on drying over silica gel in vacuum	max.1.0%
* Assay (HPLC)	min.99.0%

*:Additional test performed by Wako