

SPECIFICATION

Sodium L-Tartrate

Japanese Pharmaceutical Excipients

REQUIREMENT	SPECIFICATION
* Description Identification Optical rotation [α]D ₂₀ pH Purity (1) Clarity and color of solution (2) Sulfate (3) Oxalate (4) Heavy metals (5) Arsenic Loss on drying Assay (as anhydrous) (after drying)	Colorless crystals or white crystalline powder to pass test +25.0~+27.5° 7.0~9.0 - to pass test not more than 0.019% to pass test not more than 10ppm not more than 2ppm 14.0~17.0% not less than 98.5%

Japanese Pharmaceutical Excipients(*:Reference item)