

Code No.19-1617

# SPECIFICATION

## Saikosaponin b2

for the Japanese Pharmacopoeia Crude Drugs Test (for Assay and TLC)

REQUIREMENT	SPECIFICATION
Description (JP method)	White crystals or crystalline powder
Identification (JP method)	to pass test
Purity (JP method)	-
Related substances (JP method)	to pass test
Unity of peak (JP method)	to pass test
Assay (qNMR) (JP method)	under investigation
* Solubility in methanol	to pass test
* Loss on drying over silica gel in vacuum	max.10.0%

This product is suited for Saikosaponin b2 for assay and thin-layer chromatography of general tests on The Japanese Pharmacopoeia. (\*:Additional test performed by Wako)