

## SPECIFICATION

## Sinomenine

for the Japanese Pharmacopoeia Crude Drugs Test (for Assay and TLC)

REQUIREMENT	SPECIFICATION
Description	White or pale brown, crystalline powder
Identification	to pass test
Purity	-
Related substances	to pass test
Unity of peak	to pass test
Assay (qNMR)	measured value
* Solubility in methanol	to pass test
* Loss on drying over silica gel in vacuum	max.1.0%
* Assay (HPLC)	min.98.0%

This product is suited for Sinomenine for assay (for assay 2) and thin-layer chromatography of general tests on The Japanese Pharmacopoeia. (\*:Additional test performed by Wako)