

SPECIFICATION

L(+)-Tartaric Acid

Wako 1st Grade

REQUIREMENT	SPECIFICATION
Appearance	White crystals
Solubility in water	to pass test
Solubility in ethanol	to pass test
Specific rotation [α]D ₂₀ (c=20,H ₂ O) (after drying)	+12.0~+12.5°
Chloride (Cl)	max.0.001%
Sulfur compounds (as SO ₄)	max.0.05%
Heavy metals (as Pb)	max.0.001%
Arsenic (As)	max.1ppm
Iron (Fe)	max.0.001%
Assay	min.99.0%