

Code No.20-0145

SPECIFICATION

Thymolphthalein

Guaranteed Reagent

REQUIREMENT	SPECIFICATION
* Appearance Solubility in ethanol Absorbance (5mg/L, pH10.5) Loss on drying at 105°C Residue after ignition (as sulfate) Visual transition interval (pH)	White ~ slightly yellow, crystalline powder ~ powder to pass test min.0.35 max.1.0%(mass/mass) max.0.1%(mass/mass) 8.6~10.5(colorless ~ blue)

JIS K 8642:2021(*:Additional test performed by Wako but not required by JIS)