

Code No.20-0145

## SPECIFICATION

## Thymolphthalein

Guaranteed Reagent

REQUIREMENT	SPECIFICATION
* Appearance	White ~ slightly yellow, crystalline powder ~ powder to pass test  min.0.35  max.1.0%(mass/mass)  max.0.1%(mass/mass)  8.6~10.5(colorless ~ blue)
Solubility in ethanol	
Absorbance (5mg/L, pH10.5)	
Loss on drying at 105°C	
Residue after ignition (as sulfate)	
Visual transition interval (pH)	

JIS K 8642:2021(\*:Additional test performed by Wako but not required by JIS)