

## SPECIFICATION

Toluene  
for Chromatography

REQUIREMENT	SPECIFICATION
Appearance	Colorless clear liquid
Density (20°C)	0.864~0.868g/mL
Refractive index 20°C	1.496~1.498
Water	max.0.02%
Residue after evaporation	max.0.001%
Acidity (as HCl)	max.0.001%
Alkalinity (as NaOH)	max.0.001%
Sulfur compounds	to pass test
Thiophene	to pass test
Substances darkened by sulfuric acid	to pass test
Assay (capillary GC)	min.99.7%