

## SPECIFICATION

## 2,2,4-Trimethylpentane

for HPLC

REQUIREMENT	SPECIFICATION
Appearance	Colorless clear liquid
Density (20°C)	0.690~0.693g/mL
Refractive index 20°C	1.390~1.393
Absorbance (210nm)	max.0.50
Absorbance (220nm)	max.0.20
Absorbance (230nm)	max.0.05
Absorbance (254~400nm)	max.0.01
Water	max.0.01%
Residue after evaporation	max.0.001%
Acidity (as CH <sub>3</sub> COOH)	max.0.001%
Peroxides (as H <sub>2</sub> O <sub>2</sub> )	max.5ppm
Fluorescence test	to pass test
Assay (capillary GC)	min.99.7%