

SPECIFICATION

0.5mol/L Tetrabutylammonium Hydroxide Solution
for High Performance Liquid Chromatography

REQUIREMENT	SPECIFICATION
Appearance	Colorless ~ slightly yellow, clear liquid
Absorbance (240nm)	max.0.50
Absorbance (254nm)	max.0.30
Absorbance (300nm)	max.0.15
Absorbance (350nm)	max.0.10
Bromide (Br)	max.2.0%
Fluorescence test	to pass test
Concentration	0.45~0.55mol/L