

Code No.20-0874

SPECIFICATION

Tetrahydrofuran

Guaranteed Reagent

REQUIREMENT	SPECIFICATION
* Appearance Assay (C ₄ H ₈ O) (GC) Density (20°C) Refractive index 20°C Water Peroxides (as H ₂ O ₂)	Colorless clear liquid min.99.5%(area/area) 0.884~0.889g/mL 1.406~1.409 max.0.05%(mass/mass) max.0.005%(mass/mass)

JIS K 9705:2024*:Additional test performed by Wako but not required by JIS(with antioxidants)