

Code No.20-0874

SPECIFICATION

Tetrahydrofuran

Guaranteed Reagent

REQUIREMENT	SPECIFICATION
* Appearance	Colorless clear liquid
Assay (C ₄ H ₈ O) (GC)	min.99.5%(area/area)
Density (20°C)	0.884~0.889g/mL
Refractive index 20°C	1.406~1.409
Water	max.0.05%(mass/mass)
Peroxides (as H ₂ O ₂)	max.0.005%(mass/mass)

JIS K 9705:2024*:Additional test performed by Wako but not required by JIS(with antioxidants)