

Code No.20-0896

## SPECIFICATION

## 1,1,3,3-Tetramethoxypropane

Wako 1st Grade

| REQUIREMENT           | SPECIFICATION                                |
|-----------------------|----------------------------------------------|
| Appearance            | Colorless ~ slightly yellow,<br>clear liquid |
| Density (20°C)        | 0.991~0.998g/ml                              |
| Refractive index 20°C | 1.406~1.410                                  |
| Water                 | max.0.2%                                     |
| Assay (capillary GC)  | min.98.0%                                    |