

## SPECIFICATION

0.5mol/L Tetrabutylammonium Dihydrogenphosphate Solution  
for HPLC

| REQUIREMENT        | SPECIFICATION                             |
|--------------------|---|
| Appearance         | Colorless ~ slightly yellow, clear liquid |
| Absorbance (240nm) | max.0.50                                  |
| Absorbance (254nm) | max.0.30                                  |
| Absorbance (300nm) | max.0.15                                  |
| Absorbance (350nm) | max.0.10                                  |
| Bromide (Br)       | max.0.1%                                  |
| Fluorescence test  | to pass test                              |
| Concentration      | 0.45~0.55mol/L                            |