

SPECIFICATION

2mol/L Triethylamine Acetate Solution, pH7.0
for Genetic Research

REQUIREMENT	SPECIFICATION
Appearance	Colorless clear liquid
Absorbance (254nm)	max.0.030
Absorbance (260nm)	max.0.015
Absorbance (270nm)	max.0.010
Absorbance (300~350nm)	max.0.010
pH	6.9~7.1
Aldehydes (as CH ₃ CHO)	max.0.0020%
Diethylamine	max.0.0050%
Ethylamine	max.0.0010%
Acetic acid	1.9~2.1 mol/L
Practical test	to pass test