

SPECIFICATION

25% Tetramethylammonium Hydroxide Solution
Super Special Grade

REQUIREMENT	SPECIFICATION
Appearance	Colorless clear liquid
Chloride (Cl)	max.0.5ppm
Sodium (Na)	max.0.01ppm
Potassium (K)	max.0.01ppm
Copper (Cu)	max.0.005ppm
Calcium (Ca)	max.0.01ppm
Zinc (Zn)	max.0.01ppm
Cadmium (Cd)	max.0.003ppm
Mercury (Hg)	max.0.003ppm
Lead (Pb)	max.0.005ppm
Arsenic (As)	max.0.01ppm
Chromium (Cr)	max.0.003ppm
Manganese (Mn)	max.0.005ppm
Iron (Fe)	max.0.005ppm
Cobalt (Co)	max.0.003ppm
Nickel (Ni)	max.0.005ppm
Assay	24.0~26.0%