

## SPECIFICATION

 $(\pm)$ - $\alpha$ -Tocopherol Acetate Standard  
for Food Analysis

REQUIREMENT	SPECIFICATION
Appearance	Colorless ~ yellow, clear liquid
Solubility in ethanol	to pass test
Refractive index 20°C	1.494~1.499
Specific absorbance E1%1cm (284nm)	41.0~45.0
Water	max.0.2%
$\alpha$ -Tocopherol	max.0.5%
Assay (HPLC)	min.96.0%