

SPECIFICATION

Tetrahydrofuran, Deoxidized, Stabilizer Free
for Organic Synthesis

REQUIREMENT	SPECIFICATION
Appearance	Colorless clear liquid
Density (20°C)	0.884~0.889g/mL
Refractive index 20°C	1.406~1.409
Absorbance (250nm)	max.0.50
Absorbance (280nm)	max.0.20
Absorbance (350nm)	max.0.02
Water	max.0.001%
Dissolved oxygen	max.1ppm
Peroxides (as H ₂ O ₂)	max.0.005%
Assay (capillary GC)	min.99.5%