

Code No.20-1915

## SPECIFICATION

Tilmicosin Standard(mixture of isomers)  
for High Performance Liquid Chromatography

| REQUIREMENT  | SPECIFICATION                              |
|--|--|
| Appearance   | White ~ nearly white, crystalline powder ~ |
| Solubility in methanol   | powder                                     |
| Specific rotation [ $\alpha$ ] <sub>D20</sub> (c=0.1,CH <sub>3</sub> OH) | to pass test                               |
| Water  | under investigation                        |
| Assay (HPLC) (mixture of isomers)  | under investigation                        |
|  | min.98.0%                                  |