

## SPECIFICATION

Tetrahydrofuran, with Stabilizer  
for GPC

REQUIREMENT	SPECIFICATION
Appearance	Colorless clear liquid
Density (20°C)	0.884~0.889g/ml
Refractive index 20°C	1.406~1.409
Absorbance (300nm)	max.1.00
Absorbance (320~600nm)	max.0.10
Water	max.0.05%
Peroxides (as H <sub>2</sub> O <sub>2</sub> )	max.0.005%
Suitability for GPC test	to pass test
Assay (capillary GC)	min.99.5%
Stabilizer (BHT)	0.01~0.03%