

SPECIFICATION

Fluorescein Sodium

Guaranteed Reagent

REQUIREMENT	SPECIFICATION
* Appearance	Green or yellowish red ~ reddish brown, crystalline powder ~ powder or mass to pass test
Solubility in water	
Ethanol-insoluble matter	max.1.0%(mass/mass)
Absorbance (3mg/L, pH7.0)	min.0.45
Loss on drying at 105°C	max.7.0%(mass/mass)
Residue after ignition (as sulfate) (after drying)	35.0~40.5%(mass/mass)
Chloride (Cl)	max.0.1%(mass/mass)
Zinc (Zn)	max.0.05%(mass/mass)
Sensitivity	to pass test

JIS K 8830 : 2023(*:Additional test performed by Wako but not required by JIS)