

Code No.22-0037

SPECIFICATION

Vitamin K1

REQUIREMENT	SPECIFICATION
Appearance	Yellow ~ reddish yellow, clear liquid
Solubility in 2,2,4-trimethylpentane	to pass test
Ratio of absorbance(E253.5nm/E248.5nm)	0.70~0.75
Ratio of absorbance(E253.5nm/E269.5nm)	0.74~0.79
Heavy metals (as Pb)	max.0.002%
Menadione	to pass test
Assay	97.0~102.0%