

## SPECIFICATION

Wakogel(R) C-100  
for Column Chromatography

REQUIREMENT	SPECIFICATION
Appearance	White small granules
Particle size (150~425 $\mu\text{m}$ )	min.75%
Particle size (coarser than 425 $\mu\text{m}$ )	max.3%
Precipitation volume	1.8~2.4ml/g
Loss on drying at 180°C	max.5.0%
pH (50g/l aqueous slurry,25°C)	5.5~7.0
Flowing rate	min.3.5ml/min
Separation power	to pass test
Activity	to pass test