

## SPECIFICATION

Wakogel(R) Q-50, 75~250  $\mu$ m(60~200mesh)  
for Column Chromatography

REQUIREMENT	SPECIFICATION
Appearance	White powder
Particle size (75~250 $\mu$ m)	min.90%
Precipitation volume	1.2~1.7ml/g
Loss on drying at 150°C	max.5.0%
pH (50g/l aqueous slurry,25°C)	4.0~6.0
Practical test (Separation power)	to pass test
Practical test (Polymerization power)	to pass test