

SPECIFICATION

Wakogel(R) S-1
for Column Chromatography

| REQUIREMENT | SPECIFICATION |
|-----------------------------------|---------------|
| Appearance | White powder |
| Particle size (63~425 μ m) | min.70% |
| Precipitation volume | 3.2~5.8ml/g |
| Loss on drying at 130°C | max.5.0% |
| pH (20g/l aqueous slurry,25°C) | 4.5~6.5 |
| Hexane elution test | to pass test |
| Suitability for PCB determination | to pass test |
| Flowing rate | min.1.4ml/min |