

Code No.23-0143

## SPECIFICATION

Wakogel(R) C-200E  
for Column Chromatography

REQUIREMENT	SPECIFICATION
Appearance	White, crystalline powder or small granules
Particle size (75~150 $\mu$ m)	min.75%
Particle size (pass 45 $\mu$ m)	max.5%
Precipitation volume	1.8~2.7ml/g
Loss on drying at 180°C	max.6.0%
pH (50g/l aqueous slurry,25°C)	5.5~7.5
Flowing rate	1.6~3.0ml/min
Separation power	to pass test
Activity	to pass test