

## SPECIFICATION

WakosilR C-300  
for Column Chromatography

REQUIREMENT	SPECIFICATION
Appearance	White powder
Particle size (38~63 $\mu\text{m}$ )	min.80%
Particle size (coarser than 63 $\mu\text{m}$ )	max.5%
Particle size (finer than 38 $\mu\text{m}$ )	max.10%
Precipitation volume	1.7~2.3mL/g
Loss on drying at 180°C	5.0~10.0%
pH(50g/L aqueous slurry, 25°C)	6.5~7.5
Flowing rate	1.2~1.8mL/min
Separation power	to pass test
Activity	to pass test