

Code No.23-0232

# SPECIFICATION

## Wogonin

for the Japanese Pharmacopoeia Crude Drugs Test (for Thin-Layer Chromatography)

REQUIREMENT	SPECIFICATION
Description (JP method)	Yellow crystals or crystalline powder
Identification (JP method)	to pass test
Purity (JP method)	-
Related substances (JP method)	to pass test
* Loss on drying over silica gel in vacuum	max.1.0%
* Assay (HPLC)	min.98.0%

\*:Additional test performed by Wako